





Completing bacterial genome assemblies with multiplex MinION sequencing

Ryan R Wick et al. 2017



Sample	Illumina-only assembly graph	ONT reads	Hybrid assembly graph
INF042			

```
refdata=/home/abaryiames/Lab_Projects/polishtest/GCF_000240185.1_ASM24018v2_genomic.fna
reads=/usr/share/data/proj_data/pneumonia/5170843/barcode01.fastq
i_1=/usr/share/data/proj_data/pneumonia/5170831/barcode01_1.fastq.gz
i_2=/usr/share/data/proj_data/pneumonia/5170831/barcode01_2.fastq.gz
/all_segments.fasta
sample_name="barcode01"
mkdir $sample_name
cd $sample_name
```

```
logfile='log.txt'
touch $logfile
echo $datapath > $logfile
```

```
ls -lh $datapathO
```

```
#Porechop - Gets rid of adapter seqs
porechop -i $datapathO/$sample_name* -o barcode01_chop.fastq
```

```
chopped_pore=/home/abaryiames/Lab_Projects/pneumonia_script/barcode01_chop.fastq
```

```
#Canu - Assembly
canu -p canu -d canu_dir genomeSize=5.5m maxThreads=8 -nanopore-raw $chopped_pore
oxford_canu=/home/abaryiames/Lab_Projects/pneumonia_script/barcode01/canu_dir/canu.contigs.fasta
```

```
#Minimap / Racon - Polishing
minimap -Sw5 -L100 -m0 -t2 $assembly $reads > mapIT01.paf
racon -t 16 $reads mapIT01.paf $oxford_canu > mapIT01.fa
minimap -Sw5 -L100 -m0 -t2 mapIT01.fa $reads > mapIT02.paf
racon -t 16 $reads mapIT02.paf mapIT01.fa > mapIT02.fa
minimap -Sw5 -L100 -m0 -t2 mapIT02.fa $reads > mapIT03.paf
racon -t 16 $reads mapIT03.paf mapIT02.fa > mapIT03.fa
```

```
#Hybrid Assembly
unicycler -1 $i_1 -2 $i_2 -l $chopped_pore -o unicycler_hybrid_assemblies --threads 16
```

```
#Illumina Assembly
unicycler -1 $i_1 -2 $i_2 -o unicycler_illumina_assemblies --threads 16
```


```
#Name assemblies
assemblyO=/home/abaryiames/Lab_Projects/polishtest/mapIT03.fa
assemblyI=/home/abaryiames/Lab_Projects/pneumonia_script/unicycler_illumina_assemblies/assembly.f
asta

assemblyH=/home/abaryiames/Lab_Projects/pneumonia_script/unicycler_hybrid_assemblies/read_align
ment/all_segments.fasta
```


```
#dnadiff
dnadiff -p oxford_pilon $refdata $pilon
dnadiff -p oxford_pilon $refdata $assemblyI
dnadiff -p oxford_pilon $refdata $assemblyH
```


```
#Quast.py
quast.py $assemblyO -R $refdata --nanopore $reads -o quast_oxford --threads 16
quast.py $assemblyI -R $refdata --pe1 $i_1 --pe2 $i_2 -o quast_illumina --threads 16
quast.py $assemblyH -R $refdata --pe1 $i_1 --pe2 $i_2 --nanopore $reads -o quast_hybrid --threads 16
```


Read lengths vs Average read quality plot





Depth of Assemblies














1.00x

0.733x

0.613x